

Test changes

Test name: Aldosterone

Test #: 22/82088.1

Notification date: January 17, 2019

Effective date: January 22, 2019

EXPLANATION

Mayo Clinic Laboratories (MCL) has announced a change to the preferred transport temperature for the Aldosterone (22/82088.1), effective January 22, 2019

DETAILS

On January 22nd, the preferred temperature for Aldosterone specimens will be updated as indicated below. The Allina Health Laboratory test catalog will be updated to reflect the changes at that time.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	30 days
	Ambient	4 days
	Refrigerated	28 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Ambient	4 days
	Frozen	30 days